Choose your Configuration

RETImap[®] Model 520

Modules		Programs
SLO Basis Unit	•	Pattern VEP
ERG-Module	О	Pattern ERG
OCT-Module	0	Pattern VEP
FA -Module	0	Flash ERG
	<u> </u>	Focal ERG
Stimulators		EOG
Stimulator Projector (DLP)	•	Multifocal E
Ganzfeld Q450	О	Multifocal \
Ganzfeld Q450SC	0	Multifocal s
Ganzfeld 0450SCX	0	Glaucoma S
		Visual Acuit
19" Stimulator Monitor	O	Nystagmog
BABYflash	О	ON-OFF Res
MINIganzfeld	О	Photopic ne
Pattern Handheld	О	S-Cone ERG
Ameril: Generate		Pupillometr
2 Channels		Calibratio
4 Channels	0	Automatic G
6 Channels	0	Automatic St

Part No	• Standard O Option
Automatic Stimulator Monitor calibration	0
Automatic Ganzfeld calibration	О
Calibration tool	
Pupillometry	О
S-Cone ERG	0
Photopic negative response	О
ON-OFF Response	О
Nystagmography	О
Visual Acuity	0
Glaucoma Screening	Ο
Multifocal science	О
Multifocal VEP	О
Multifocal ERG	О
EOG	О
Focal ERG	О
Flash ERG	О
Pattern VEP	О
Pattern ERG	О
Pattern VEP	O

Standard Accessories		Part No
	DTL Electrode S (sterile) and connection cable Pack of 2 electrodes Pack of 2 cable	Part No 1000-510-336 Part No 1000-510-304-D
777	EEG Gold-Electrodes Ø 10 mm Electrode, DIN connector 1,5mm Set 2 x red, 2 x blue, 1 x black DIN connector 1,5 mm	Part No 1000-316-300-D Part No 1000-316-301-D
Ten20	TEN20 Conductive and Adhesive Paste 230 g Paste 115 g Paste	Part No 1000-361-300 Part No 1000-362-300
uj.Nuprep 1	NuPrep Gel for Skin Preparation 115 g Gel	Part No 1000-370-300

Ο

Distributor:

8 Channels



Friedrich-Franz-Str.19, D-14770 Brandenburg, Germany Phone: +49 3381 890-1034, Fax: +49 3381 890-2994 www.roland-consult.de, e-mail: info@roland-consult.de

17-51-1.3E

Certification: Quality Management System EN ISO 13485 **TÜV Rheinland**

Made in Germany

SLO



R **RETI**map Model 520

Electrophysiology and Imaging

ERG







ALL IN ONE

Dark-adapted 0.01 ERG b-wave: rod-initiated on pathways

Dark-adapted 3.0 ERG a-wave: photoreceptor & postreceptoral on pathways b-wave: on & off bipolar cells

Dark-adapted 10.0 ERG a-wave: photoreceptor & postreceptoral on pathways b-wave: predominantly rod bipo lar cells (on pathways)

Dark-adapted 3.0 oscillatory potentials On & off pathways reflecting middle retinal layers & vascula function

Light-adapted 3.0 ERG a-wave: cones with post-receptoral on & off pathways b-wave: on & off bipolar cells

Light-adapted 3.0 fli Cone systems with ost-recentoral or







OCT

Product overview



Stimulator Monitor



Ganzfeld

Tilting Stand



Ganzfeld ERG with different colours



Features:

- Delivered with standard ISCEV programs
- Additional programs Visual Acuity, Glaucoma Screening, S-Cone ERG, Photopic negative response
- Possibility to create own programs
- Automated impedance test can start every time
- Automatically analyze with placing the markers during the examination
- Automatically artifact elimination in all programs, also for mfERG and mfVEP
- Ganzfeld with 5 different colors and eye fixation camera
- Possibility to measure the pupil size automatically
- Optimized short screening tests for children
- Digital Filter for signal processing
- Show typical curve in the analyze and printout mode
- Delivery with normal values, it is easy to insert the own normal values
- Export all data to EXCEL
- Working in the LAN, Reader Stations
- Service via Team Viewer

Operating unit:

- Taking account of the existing technological level
- Intel core i5-2500 CPU (3.3 GHz), 2 GB RAM, HDD: 500 GB
- 24"TFT control monitor, keyboard
- Mouse, colour printer
- Table dimensions (W x H x D) $120 \times 95\pm 20 \times 70$ cm
- Software : Roland Software based on MS Windows, Nero, Antivirus, Team Viewer

Bio signal amplifier:

- 2, 4, 6 or 8 channel
- Type BF, check voltage 1,5 kV
- Impedance 2 x 100 M Ω , noise < 4 μ V(SS)
- Common mode rejection >120 dB
- Gain up to 150.000
- Sensitivity 10µV/div to 2 mV/div
- High pass: 20 Hz to 10 kHz
- Low pass : 0,02 Hz to 1 kHz

Stimulator monitor

- 19"-color monitor, max. luminance 220 cd/m²;
- high contrast

MINIganzfeld

- Checkerboards, bars fields: full, half or quarter
- Pattern reversal / appearance/disappearance
- Software controlled contrast settings (3 % 99 %)
- black and white or different colour settings
- variable fixation points,
- special pictures for children









- The RETImap Model 520 is a new high technology combination of electrophysiology and fundus imaging with OCT and Angiography. It's a modular system and the basic unit consists of:
- infrared confocal Scanning Laser Ophthalmoscope cSLO
- integrated DLP stimulator unit

and can be combined with:

- Electroretinogram (ERG) Module + Visual Evoked Potentials (VEP) Module
- Fluorescein Angiography (FA) Module +Autofluorescence (AF) Module
- Spectral Domain Optical Coherence Tomography (SD-OCT) Module
- Amplifier 2,4,6, or 8 channels
- Flash Stimulators: Ganzfeld, BABYflash, MINIganzfeld, Pattern Handheld
- Pattern Stimulators: 19" TFT Monitor, Pattern Handheld (LED)

The basic unit delivers a highly detailed image of the fundus. The integrated Roland Consult Fundus Eye Tracker (RCFET) supports to find exactly the same position of past examinations. The fixation target is shown by a high resolution DLP projector. Colours and design of the target can be customizable. Possibility of standard programs:

- Optical Coherence Tomography (OCT)
- Fluorescein Angiography, Auto fluorescence
- focal ERG, mf ERG
- focal VEP, mf VEP
- Pattern ERG,
- Flash ERG (Ganzfeld 6 steps),
- ON-OFF Response, Photopic Negative Response, S-Cone ERG
- EOG (fast, slow)
- Visual Acuity
- Glaucoma Screening Contrast Sensitivity
- Nystagmography
- Pupillometry
- More different examinations are freely customizable.

The software includes many advanced features like FFT analyse, OFF-line-averaging and pre-trigger. All electrophysiological results can exactly overlay to the SLO image. Different combination of structure and function results in one protocol.

Ganzfeld Q450

The Ganzfeld consists of a 400 mm full field globe, with a central fixation LED and two EOG fixation LEDs. The brightness of these LEDs are software controlled and an infrared camera is integrated. There are two models Q450C and Q450SC. The Xenon tube module (X) for high intensity flash is optional available for both models. Model Q450C(X): white, blue, red Model Q450SC(X): white, blue, red, royal blue, green, amber Stimulus ON-OFF Stimulus white flash - standard intensity 3,0 cds/m² - Interval -40 dB to +5 dB in steps of 5 dB Stimulus colour flash: - standard flash 3,0 cds/m² - royal blue (455 nm) interval -50 dB to -5 dB in steps of 5 dB – blue (470 nm) interval -45 dB to 0 dB in steps of 5 dB (525 nm) interval -45 dB to 0 dB in steps of 5 dB - green - amber (590 nm) interval -45 dB to 0 dB in steps of 5 dB (625 nm) interval -45 dB to 0 dB in steps of 5 dB – red Stimulus ON-OFF: - all colours: 1 ms to 1000 ms adjustable in steps of 1 ms Background Luminance: adjustable in 1,0 cd/m² steps 1000 cd/m² - white: - royal blue (455 nm): 100 cd/m² (470 nm): 200 cd/m² - blue – red (625 nm): 200 cd/m² (525 nm): 500 cd/m² - green - amber (590 nm): 750 cd/m² simultaneous LED use to generate different flash/background colours Option X: - Xenon tube for white high flash - Intensities: 9,5 cds/m² (+5 dB), 30 cds/m² (+10 dB), 95 cds/m² (+15 dB) For model Q450 SC(X) are additional flexible settings possible. Option Flimmer check: (according Prof. Kremers) Only for model Q450 SC(X) For each colour: - Selectable waveform type: sine wave, rectangular wave, - Triangular wave, ramp up or ramp down - Phase shift: 0°- 359° in steps of 1° - Contrast: 0,1 % -100 % in steps of 0,1 % - Stimulation frequency: 1 Hz - 150 Hz Option: Pupillometry - Full field Ganzfeld stimulation - Resolution time 33 ms (30 images per second) - Resolution pupil size 0.1 mm - Examination settings number of cycles, cycle time, record time, flash time, flash intensity, averaging of cycles **OPTIONS** MINIganzfeld - Flash Luminance: standard flash: 3,0 cds/m² white, - Range -25 dB...+5 dB in steps of 5 dB - Backlight: 10, 30, 50, 100, 450 cd/m² Pattern Handheld - Checkerboard: 10 x10 mm, fields á 5 x 5 mm - Luminance 80 cd/m²; Pattern Reversal BABYflash - Standard flash: 3,0 cds/m² white, blue, red: range -25 dB...+10 dB - Backlight: 10, 30, 50, 100, 450 cd/m² Calibration Tool - Measurement with Mayo Monitor via USB connection - Automatic Ganzfeld calibration Automatic Stimulator monitor calibration

Moduls

cSLO image with OCT line scan

SLO

Scanning Laser Ophthalmoscope Unit (SLO) It allows basic examinations and has following

Field of view: $30^{\circ} \times 30$ Laser source:SLD 920Optical resolution: $<30 \ \mu m$ Focus adjustment: $-15 \ to +$ Digital image: $12 \ bit$ Image Resolution: 512×51 Pupil size requirement: $\geq 2 \ mm$ EyeTracker:software

30° x 30° SLD 920 nm (IR Image) <30 µm -15 to +15 diopter (D) 12 bit 512 x 512 :≥ 2 mm

software module

ERG VEP

Module Electroretinogram (ERG) and Visual Evoked Potentials (VEP)

It allows simultaneous infrared laser monitoring during electrophysiological diagnostic. With the built-in stimulator, the high resolution DLP projector, it is possible to stimulate a large number of retinal locations and extract their responses.

The final result is a function map superimposed onto the fundus cSLO image. Regionally confined areas of dysfunction can be detected. The Roland Consult Fundus EyeTracker (RCFET) software helps to detect all eye movement artefacts during the test. This device can also be upgraded with the classical electrophysiological stimulators and programs. The various electrophysiological functions are:

- focal ERG
- multifocal ERG
- Pattern ERG
- Flash VEP
- Pattern VEP
- multifocal VEP

Field of view:	25° x 25°
Monitor stimulator:	DLP projector
Brightness:	300 cd/m ²
Vertical frequency:	60 Hz
Resolution:	800 x 600 Pixel
Color Stimulator RGB:	636nm/520nm/452nm
Stimulator Protocols:	focal ERG, multifocal ERG,
	Flash VEP and Pattern VEP
Control monitor:	10 inches



FA

Module Fluorescin Angiography (FA)

The RETImap Angiography module runs with a blue laser source 488 nm. It generates higher quality fundus images than alternative white light illumination photo systems, and since no flashes are involved, the system is more patient friendly. At the beginning of an injection the laser source automatically switches to Angiography mode and a video is recorded with an exact timestamp on every image. This helps the user to concentrate more on the patient during the examination and review the angiography images later. Another advantage of the laser system is the high speed of image acquisition. Instead of watching static images from early, mid, or late phase, it is easy to observe the dye in the vessels to localize narrowings and partial blockades in real time.

Laser source Blue:488 nmDigital image:512 x 512Record Mode:15 pictures / sec

OCT

Module Optical Coherence Tomography (OCT)

The combination of RETImap cSLO + OCT enables the exact analysis of the morphology. The user just selects an area on the cSLO image and describe the scan (hor./ ver. LineScan, CircleScan, VolumeScan). Because of an optimized hardware design a single OCT Scan originally comes in a very high quality. In combination with the Roland Consult Robust OCT Tracker (RCROT) the signal noise ratio will be greatly reduced to a minimum. The OCT software is able to detect up to 7 layers automatically and visualize them in thickness maps. The integrated DLP Beamer helps the patient with an internal customizable fixation cross and reduces the motion artefacts.

Laser Source: Max. scan speed: Scan depth: Axiale resolution: Transverse Resolution:

830±50 nm, Super Luminescent Diode (SLD) 25 000 scan / sec. 2-2,5 mm 5 μm on: 15 μm



cSLO image with OCT LineScan



Angiography video with different phases









mfERG results superimpose on cSLO image



automatic layer detection





Measurement and Results



mfERG

mfVEP

N75 N135 N75 🖌 ¥N135

normal mfERG plots

pattern VEP



mfVEP plots





Pattern Handheld



BABYflash ERG/VEP



MINIganzfeld for flash ERG/VEP











Ganzfeld Q450 EOG stimulus







EOG result



analyse regression curve



EOG printout normal result



pupillometry





nystagmography measurement left



steady state VEP



steady state photopic 30 Hz ERG









overlapping with visualfield

र ज म म

pattern ERG

3D plot normal patient

Contra Co

w









transient VEP



transient VEP



photopic ERG



pupillometry measurement

pupillometry result



nystagmography results

