



Xltek® Ambulatory EEG/PSG

Patient Comfort, Value, Performance & Flexibility

Comfort: Lightweight and compact amplifier comfortable for adults and children

Multipurpose: Support for EEG and PSG video studies in both clinical and ambulatory settings

Standard components: Standard off-the-shelf AA batteries and high-speed USB 2.0 connectivity

Performance: Up to 96 hours of ambulatory recording; less than 5 minutes to upload 24 hours of EEG data

Innovation: Unique ambulatory video technology enables up to 72 hours of wirelessly synchronized HD video with simplified setup in patient's home

Tools: Microsoft® SQL Server database and integrated Spike & Event detection algorithm for maximum efficiency




Now featuring
HD Video

Convenience, patient comfort & reliability for a broad range of EEG & PSG recording needs

High-Definition
Video

QEEG Trending



Small Amplifier — Big Performance

While the Xltek Trex HD device is small, it's big on the performance, flexibility and reliability you demand. It features 24 referential, 4 differential, 4 DC and Oximetry channels, enabling use for a variety of EEG and PSG video studies in ambulatory or clinical mode. Connecting to an acquisition computer through a standard high-speed USB 2.0 cable eliminates the need for expensive proprietary cables and allows uploading of 24 hours of EEG data from the amplifier to your station in less than 5 minutes.

Reliable Video EEG at a Patient's Home — now in HD video

Frequently, the best way to capture a sleep, neurological or cardiac event is to use an ambulatory device in a patient's home environment. But will the device be durable enough? Can you depend on its signal quality? Will the battery last? Is it comfortable enough for a child to wear? Xltek Trex HD, the ambulatory video amplifier from Natus, answers "Yes" to all these questions.


The Trex HD amplifier is among the lightest in the industry, weighing 2/3 of a pound and utilizing standard cost-effective AA batteries. It integrates seamlessly with Natus® NeuroWorks®/SleepWorks™ software to provide fully synchronized EEG and video data. It is easy to set up in a patient's home, requiring no cabling and no need for the patient to interact with the software, and is capable of recording up to 96 hours of video and audio data without the need to free storage space. The patient can start and stop video while continuously recording EEG and without affecting video synchronization. This provides you with the convenience and cost-effectiveness of home ambulatory monitoring complemented with high quality video playback.





Amplifier features

- Compact (4 x 6 x 1 in) and lightweight (0.66 lb)
- Ultra low noise with superb signal quality
- Ample referential, differential and DC channels for all your needs
- Integrated Oximetry channels
- Wireless synchronization with HD video
- Fast upload of ambulatory data allowing quick turnaround
- Battery indicator LEDs
- Automatic study recording after battery replacement
- Synchronized patient event marker



System features

- Desktop or laptop-based acquisition stations
- Microsoft SQL® Server database enabling quick and efficient data retrieval
- Integrated Spike & Event detection and trending for efficient data review
- Capability of running EEG and PSG tethered or untethered studies
- Synchronized HD video recording and playback
- Simple annotation entry
- Easy archiving to a network device or any type of recordable media
- Photic stimulation integrated within software



Data Security

Set up individual user or group access and permissions for the system and within the application. Assign separate lock-down rules for different computers on the same network. Access Audit Logs to review user activity on the system. Security that is integrated with Microsoft Active Directory and peace of mind with available system encryption.



Service

Excellence in customer and technical service is a core value of Natus Neurology.

Essential Support Elements:

- Accessible and effective Technical Support
- Definitive technical documentation and knowledgeable installation teams
- Excellent replacement unit or spare part availability
- Extended warranty and service coverage programs
- Comprehensive, flexible customer training courses

Supplies

Natus Neurology offers a full range of neurodiagnostic supplies, featuring:

- Dedicated and knowledgeable customer service
- Streamlined order processing to save you time
- Convenient online portal (US Customers only)

Natus Neuro Store – www.natusneurostore.com

To learn more about Natus Neurology Service Programs or our full line of Supplies, contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**



Natus Neurology Incorporated

3150 Pleasant View Road
Middleton, WI 53562 USA
Tel: 1-800-356-0007
1-608-829-8500
Fax: 1-608-829-8709
www.natus.com

Natus Medical Incorporated

Excel-Tech Ltd. (XLTEK)
2568 Bristol Circle, Oakville
ON L6H 5S1 CANADA
1-800-303-0306
1-905-829-5300

©2017 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.